für Handelsnamen: R 100.xx bis RK3-xxx EMDN-Code:S010401
Single Registration Number (SRN)DE-MF-000005592
Basic UDI: 4063058000000323G

EN

STERILIZATION CONTAINER SYSTEM **USER MANUAL**

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EN

A-FOREWORD

This user manual is a guiding feature for the longer-lasting use of our products and for the satisfaction of our customers. The recommendations have been prepared in accordance with the instructions from European standards. REDA Instrumente GmbH is a company that always cares about customer satisfaction and aims to produce high quality products. For this reason, while preparing the guideline, it is aimed to inform our customers about the questions and problems that are relevant to the subject.

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WARNING

Before using your product, please read this guide from start to finish and keep it accessible to users.

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CHAPTER-I
STERILIZATION CONTAINER MODELS

EN

KAPITEL-I
STERILISATIONSCONTAINERMODELLE

DE

REDA, since 1990, has advanced to the highest level to ensure the reliable use of medical devices and has become one of the world's leading producers of Sterilization container systems. REDA Sterilization Container systems have been produced in accordance with European standards. REDA Sterilization Container Systems are used to transport, store, sterilize and maintain sterile medical equipment for 6 months. Sterilization Container systems has the "CE" marking. Validation has been provided for Sterilization and stabilization of surgical instruments, implants and surgical motor sets of sterilized Container Systems which have been produced by fulfilling the requirements of the EC regulations. By producing Sterilization Container Systems in the Medical Industry, we offer quality that protects patient health. We serve our boxes and lids with 6 basic colour options all over the world which are produced with our special alloy, high quality 2 mm aluminium raw material. REDA container systems offer a greater alternative to the coiling method and also provides great return of investments. REDA quality that considers the needs of the market and its users, is preferred by many users all over the world due to its reasonable prices and product range for every user.

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Standard model container systems are designed for providing perfect steam flow and penetration for suitable loads (See 5.2) according to the standards and can be used with single use paper filter or reusable PTFE Teflon filters.

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
Figure /Abbildung - 2

1.2 Window Model

Window model Sterilization container system is named after transparent window which is produced from PPSU material and resistant to high temperature, scratches, deformation and chemicals. This window allows users to control the Sterilization container without opening the lid.

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
Figure /Abbildung - 4

NOTE

There is no hole in the lid of the window model container system to provide airflow. Therefore, window model boxes must be perforated. Filter, valve or Bio-stop system can be used in the box.

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DE

1.3 Bio-Barrier Model

Bio-Barrier model sterilisation container systems are designed with mechanical valves that do not require the use of any single-use filter in both the lid and the box.

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View of valve mechanism before Sterilization process. (See. Figure-5) (Rest status,+/- no pressure and valves closed.)

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During sterilization, the valves are opened by the effect of pre-vacuum. Thus, dry air comes out of the container. (See. Figure-6)

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During sterilization, the valves are opened in the opposite direction by the pressure effect.

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View of the introduction of hot steam into the container for sterilization. (See. Figure-7)

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The re-opening of the valves with va-cuum effect during the sterilisation process. The view of hot steam flowing out of the container.(See. Figure-8)

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Closure of the valve mechanism after the completion of Sterilization. (Resting position) (See. Figure-9)

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1.4 B3 Barrier Model

1.4 B3 Barriermodell

B3 model Sterilization container systems are designed to provide different filtering mechanisms in a single sterilization container lid according to user preference. B3 model sterilization container systems can be used with 3 different filtering mechanisms. (See. Figure-10)

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

Figure / Abbildung - 10
Mounting View / Montageansicht

1.4.1 B³ Filter System (PTFE)

B3 Model filtering mechanism is used with PTFE Teflon Filter which can be changed by the user easily like in standard model filter-based Sterilization containers. The filter, which is placed in a cassette, used in this system is provided only by REDA.(See. Figure-11)

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
Figure /Abbildung - 11

1.4.2 B³ Valve System

B3 Model valve system is a modernized new model of Bio- Barrier systems and has the same working principal.(See. Article 2.3) This system is based on the principle of operation of valves with the effect of pressure and vacuum during sterilization in the autoclave. (See. Figure-12)

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
Figure /Abbildung - 12

1.4.3 B3 Bio-Stop System

B3 Model Bio-Stop system is a flow inhibition system based on a principle discovered by Louis Pasteur. With a separation efficiency of 99.99 %, it offers significantly better protection against microbes and particles. (See. Figure-13)

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


Figure / Abbildung - 13

1.5 Standard Plus Model

Standard Plus model container systems are designed for providing perfect steam flow and penetration for suitable loads (See 5.2) according to the standards and can be used with single use paper filter or reusable PTFE Teflon filters.

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1.6 Plasma Model

This model is manufactured in accordance with H₂O₂ sterilization. The surface is brighter than standard container models. In this model, Box and Lid are designed as Perforated model in order to make rays easier to pass through.

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1.7 A1 and A2 Model

REDA A1 and A2 model sterilization container systems basically have similar features with standard model sterilization container systems. It has been designed by restricting some features of standard model which will not risk any of sterilization necessities.

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CHAPTER-II
COMPONENTS OF STERILIZATION CONTAINER

EN

KAPITEL-II
STERILISATIONSCONTAINERCOMPONENTEN

DE

**2.0 STERILIZATION CONTAINER
COMPONENTS AND THEIR
FEATURES**

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

Non-Perforated boxes should be used with perforated Lids for sterilization. The non-perforated box, together with the non-perforated lid, is only used for transport purposes.

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






- Grau
- Gelb
- Grün
- Blau
- Rot
- Schwarz
- Transparent

NOTE

Non-Perforated lids should be used with perforated boxes for sterilization. The non- perforated lid, together with the non-perforated box, is only used for transport purposes. (See 2.4)

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




2.4 Box, Lid and Mechanism Compliance Scheme

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

Wire Mesh-Sheet Basket
Drahtblechkorb

Perforated Sheet Basket
Perforierter Metallkorb

There are 3 kinds of basket system: Wire Basket, Sheet Metal Basket and Wire-Metal Basket. Baskets are used for the grouping of medical hand tools (set-up), washing and stacking in sterilization container.

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B

Figure /Abbildung - 21

It is enough to insert the sheet metal lid inward (in the direction of the arrow) by inserting your finger into the lock tongue as shown. You need to do the same thing to close it.(Figure 21-a) It is enough to turn the wire basket lid by inserting your finger into the lock tongue (in the direction of the arrow) as shown (Figure 21-b). You need to reverse the process to close it.

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
WARNING

It is recommended to use wire bases in perforated boxes(ex.with filter, valve, B3).

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

Wire Base
Drahtbasis Figure /Abbildung - 22

2.6 Single Use Components

The Disposable components are paper indicator labels, Disposable Paper filter, Disposable Polypropylene Filters and safety locks with indicators. (For more information see. 5.6)

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WARNING

- Disposable components should be replaced before each sterilization cycle.
- PP filter is only suitable for H₂O₂.

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2.7 Resuable Components

2.7.1 Resuable Filters- Standard Filters

Standard Filter system container models have 3 different reusable filters

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Table / Tabelle -4

2.7.2 Resuable Filters-B³ Barrier Model

B³ Barrier the filter module used in the model system has a reusable PTFE Filter

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
Figure / Abbildung - 24

WARNING

The service life of the B3 Barrier model filter module is 2000 sterilization cycles.

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2.8 Automatic Lock System

The automatic locking system, which can be used instead of the safety seal in the sterilization container systems, is mounted on the nameplate so that it is positioned next to the lid lock on the container faceplate. In the system which works with thermal sensitivity, the lock tongue is located in the inside-position before sterilization. During the sterilization, the tongue of the automatic locking system is positioned outwardly by the help of the heat applied to the container, thus preventing the opening of the container lid.

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EN

DE

The users need to push the latch inside as shown in the figure in order to open the lid after sterilization (Figure-27). This must be done when the container cools down till room temperature.

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
Figure /Abbildung - 27

WARNING

Do not unlock the automatic lock before the container reaches room temperature after sterilization. Otherwise, the automatic lock loses its function.

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2.9 Product Shelf Life

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Product Name	Shelf Life
Container (Box)	10 YEARS
Container (Lid)	
Wire Basket	
Plate Basket	
Aluminum Basket	
Auto-Lock Mechanism	
Disposal Stamp	
Label Holder	
Filter Dispenser	
Silicone Mat	
Paper Filter	5 YEARS
Polypropylen Filter	
Teflon Filter	
Silicone Mesh Mat	
Emulating Indicator Card	
Security Seal	
Silicone Handle	

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EN

CHAPTER-III **MAINTENANCE AND CONTROLS**

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3.0 Maintenance Instructions

3.0 Wartungshinweise

For sterilization container systems, check the following points before each use.

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3.1 Lid Gaskets control

3.1 Deckeldichtungskontrolle

Before each sterilization, the sealing gaskets on the edge of the lid should be visually inspected.

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
Figure /Abbildung - 28

In the event of cracks or tears on the sealing gasket before use, the container has lost its ability to sterilize. Replace the sealing gasket.
See. 3.1.1

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
Figure /Abbildung - 29

Lid should not be used if the sealing gasket is removed. Maintenance and repair are necessary. Replace the sealing gasket.
See. 3.1.1

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EN

3.1.1 Exchange of the lid gaskets

Choose the convenient cutting dimensions for the gaskets of the lids from the table below:

Container Lid Sizes Abmessungen der Containerdeckung	Gasket Cut Size Größe der Dichtung
1/1 Size	158 cm
3/4 Size	136 cm
1/2 Size	101 cm
Extra Large (XL) Size	180 cm
Extra Long Size	185 cm
Dental Size	89 cm
1/2 Dental Size	72 cm
Large Dental Size	118 cm
Mini Model Size	77 cm
Endo Model Size	130 cm

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Table / Tabelle -6

- Remove the gasket to be replaced from the lid.
- Make sure that the surface is cleaned using alcohol so that the adhesive silicone residue does not remain.
- Make sure that the new seal is cleaned from the residues on the contact surface.
- Apply silicone adhesive to the gasket channel.
- It should be left for 12 hours for drying.
- When it is dried, the replacement is completed.

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NOTE

We recommend you to use Henkel Sista® silicone.

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3.2 Lock System Control

Check the function of the locking system. Place the lid on the box. Close the locks. Make sure the lid is tightly closed so that it does not move over the box. Otherwise, the lock function is not fully fulfilled and the seal cannot be secured. It will not be suitable for sterilization.

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3.3 Filter Holder Control

Check that the Filter Holders of the Sterilization Container Systems with Filter Mechanism are Functioning Correctly before Use.

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

Figure /Abbildung - 30

As show in Figure 30, remove the filter holder from the lid by pressing the latch on it.

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Figure /Abbildung - 31

Check the gaskets on the inside edges of the filter holder. Do not use torn or deformed gaskets on the filter holders.


Figure /Abbildung - 32 Press your finger on the button until you hear the "click" sound. Try to rotate the filter holder by hand. If the filter holder does not rotate, the mechanism is seated.

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3.4 Filter Control

3.4.1 Standard Model Filter Control

Disposable filters (Paper Filter) must be replaced after each use. Check before use if the reusable filters (Textile filter, PTFE and Nonwoven PP filter) are available. Do not use the filters as damp or wet. Reusable filters should not use existing tensions, such as tears, cracks and other damage markings. For filter usage time see. Table-4

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
3.4.2 B3 Barrier Model Filter Control

For the B3 Barrier Model filter maintenance, unscrew the module as shown in Figure-33,

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
Figure /Abbildung - 33

C

Flip the lid over. Remove the module by turning the mechanism to the left. (Figure 33-A)

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Turn the module over. (Figure 33-B)

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Open the rear lid by turning the latch in the direction of the arrow as shown in Fig. 33-C)

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
Figure /Abbildung - 34

Check the filters before each usage. Do not use the filters as damp or wet. The filters should not have any tensions, tears, cracks and other damage markings. For filter usage time. (See Table-4)

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

Bio-Barrier Model
Bio-BarrierModell

B3 Barrier Model-Valve Modüle
B3 Barrierenmodell-Ventilenmodul

3.5.1 Bio-Barrier Model-Valve System Control

To check the Bio-Barrier model valve mechanism, press with a pointed material as shown in figure-35. If the mechanism is moving, it means that the system is running.

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
Figure /Abbildung - 35

3.5.2 B3 Barrier Model Valve Module Control

The control of the mechanism can be done while on the lid. Turn the lid over and press with the finger on the pin located in the middle of the mechanism. If the mechanism is moving, it is working.

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
Figure /Abbildung - 36

3.6 Box Lock latch Control

Check the lock latch in the box. The lid will not fit properly fit to boxes, where the latches are bent or broken. In this case the sealing cannot be fully achieved.

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
Figure /Abbildung - 37

WARNING

After the lid is correctly positioned on the box, the locking operation should be performed. The lid will cause unbalanced force to be applied when it is not fully seated on the box. For this reason, if the lock system fails, the seal will not function properly and there will be a risk of sterilization

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
Figure /Abbildung - 38

CHAPTER-IV
CLEANING AND DISINFECTION

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DE

4.0 Cleaning

The users of the sterilization container system should apply their preparations for reuse, the solutions to be used during the cleaning process and the way they are applied to the products should be in accordance with the instructions given by the manufacturer.

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4.1 De-assembly of components prior to cleaning


Remove all detachable parts of the container prior to cleaning.

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
- Remove paper labels, metal badges and auto lock components (if available) on the box

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For models with standard filter mechanism;

- Remove the filter holders and filters from the standard box and lids with filter.

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EN


DE

For B3 Barrier Model Filter or Valve module;

Remove the valve or filter module as shown below.

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Remove the Bio-Stop module as shown below.

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
Figure /Abbildung - 40-a

a- Turn the Bio-Stop module in the direction of the "unlock" sign. Remove the module top component.

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

Figure /Abbildung - 40-b

b- Remove the module subcomponent by turning left. Disconnect the components.

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WARNING

Do not expose the B3 model to temperatures above 90 °C when disassembled. Otherwise, the PPS components will undergo deformation. If thermal disinfection is to be carried out, first install the components in the main product.

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4.2 Washing Process

4.2.1 Pre-Cleaning

Sterilization containers should be rinsed without waiting for intensive medical tissue and blood stains on the after use. Pre-cleaning can be achieved with a cloth or soft bristle brush.

WARNING

In order to remove coarse dirt in the prewash it is absolutely not necessary to use cleaning devices which may cause scratching and hard surface, except soft bristle brushes. Scratches that occur on the surface of the container will cause the loss of antiseptic properties.

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4.2.2 Washing water quality

Tap water: It should be used for first wash (soaking) and for rinsing in process stages.

Critical water: It is obtained by removing microorganisms, inorganic and organic materials from water. It must be used in the mechanical washing process as machine water, in steam sterilization as steam water and in the last rinsing process.

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Table / Tabelle -7

If the recommended water is not used;

- Containers are coated with anodizing to obtain antiseptic properties and resistance to corrosion. Anodized surface will suffer damage over time.
- Sterilizers consume too much water during the cycle. The harder the water, the more the valves and injectors become clogged. The possibility of deterioration of the clogged injectors becomes high and the sterilization process becomes risky.

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4.2.3 Washing Process - Manuel washing

- Firstly, dismantle all detachable parts of the container as described in section 4.1.
- Wash with cold water for at least 2 minutes. (Tap water see 4.2.2). When rinsing, ensure that it is cleaned from dirt with a soft bristle brush and a lint-free cloth.
- Submerge container in Neutral pH enzymatic detergent and soak for 5 minutes at least.
- Then rinse for at least 2 minutes so that there is no solution on the surface. (Critical water see 4.2.2)

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WARNING

When rinsing, ensure that it is cleaned from dirt with a soft bristle brush and a lint-free cloth. In order to use the enzymatic and neutral cleaning solutions properly, we recommend you to follow the manufacturer's instructions for preparing the solution. (Proper Concentration, Temperature, etc.)

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EN

DE

4.2.4 Washing Process - Mechanical washing

The containers after each use, should be washed with a properly diluted, enzymatic/ Neutral pH detergent solution recommended for use on anodized aluminum.

- Firstly, dismantle all detachable parts of the container as described in section 4.1.
- Wash with cold water for at least 1 minutes. (Tap water see 4.2.2). When rinsing, ensure that it is cleaned from dirt with a soft bristle brush and a lint-free cloth.

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Table / Tabelle -8

4.2.4.1 Specific Washing Process

The parameters are specified according to UK standards.

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DE

4.2.5 Cleaning of reusable filters

- Reusable filters can be hand-washed with neutral enzymatic detergent using a soft sponge or soft bristle brush.

WARNING

Clean the filters in such a way that they do not cause damage. Reusable filters should not use existing tensions, such as tears, cracks and other damage markings. Resuable filters should thoroughly dry after cleaning. Do not use wet or damp filters.

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4.2.6 Washing Process - Points to be taken into consideration

- Do not use alkaline and acidic detergent for washing sterilization container.
- Do not use caustic (NaOH) solutions.
- Apply the detergents to be used in accordance with the instructions given by the manufacturer. (Concentration, temperature etc.)

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NOTE

In washing process applied for surgical instruments in automatic washing devices, neutralizing (acidic) detergents are used with alkaline detergent. Do not wash the containers with the surgical instrument wash program.

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5.0 Sterilization Process

5.1 Controls before sterilization

- Before sterilization, ensure that the cleaning of the sterilization container systems is complete and correct.
- Make sure that the maintenance and controls specified in Section 3 are fully performed.

WARNING

Sterilizers vary in design and performance characteristics. It is strongly recommended that the user verify the cycle parameters for the specific sterilizer prior to use.

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5.2 Material loading for sterilization

- Do not load inside the container in such a way as to cause the container holes to close when loading the materials to be sterilized.
- Keep a distance of at least 2 inches (about 5 cm) between the lid internal surface and the load.
- Because the filter mechanism found on the lid has 2,5 cm length. Minimum 1 inch(aprxt. 2,5 cm) distance should be exist between load and the filter mechanism. These recommendations also apply to loading textiles in the container.

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WARNING

If the container holes are closed, airflow during sterilization will be prevented and the container will be crushed under vacuum.

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Table / Tabelle -9

WARNING

After sterilization, pay attention to the given loading weights to avoid excessive condensation on the materials.

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Place the textile products in the container in the upright position with the fold. When the container is full, it should be in position where a flat hand can enter easily between the textile products. Load textile products so they do not lid the container holes. Otherwise, the steam will not penetrate into the container and the sterilization .

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Figure / Abbildung - 41

5.3 Loading containers into sterilizer

For steam sterilizers;

- Do not use outer packaging for the sterilization container system.
- If more than one container system is subject to sterilization at the same time, large volumes should be placed at the bottom, small volumes being at the top.

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WARNING

On the sterilizer, do not stack the fabric packaged set on top of the sterilization container.

If the container and fabric package are to be sterilized together, stack the container on the fabric packages.

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EN

DE

5.4 Suitability for steam sterilization procedures

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Table / Tabelle -10

¹ In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer table 10 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

² Only sterilization containers guaranteed evaporation or discharge of condensate are suitable above.

³ Gravitation procedures having no drying process do not guarantee of sterility a sufficient after-drying effect.

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Table / Tabelle -11

WARNING

In the sterilization of container models with valve mechanism, the working pressure of the autoclave is 2.3 bar and the vacuum time is 30 minutes.

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According to the directives of "Health Technical Memorandum 01-01: Management and decontamination of Surgical instruments (Medical Devices- U K) used in acute care Part C: Steam sterilization." parameters of sterilization process is defined as 3 minutes between the temperature range of 134-137 °C.

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5.5.2 H₂O₂ Sterilization

Sterilization validation has been carried out with EOSL VAPROX HC Sterilant.

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The sterilizer used during Validation includes the following Systems.

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- V-PRO 1/1 Plus/maX/maX 2 cycles use 2.1 g sterilant injection of VAPROX HC Sterilant (59% H₂O₂ solution) each sterilization pulse, with 4 sterilization pulses per cycle, and a 50 °C chamber temperature.
- V-PRO maX 2 Fast Non Lumen Cycle uses 2.1 g sterilant injection of VAPROX HC Sterilant (59% H₂O₂ solution) each sterilization pulse, with 2 sterilization pulses per cycle, and a 50 °C chamber temperature.
- V-PRO 60 cycles use a 1.1 g sterilant injection of VAPROX HC Sterilant (59% H₂O₂ solution) each sterilization pulse, with 4 sterilization pulses per cycle, and a 50 °C chamber temperature. The volume of the V-PRO 60 sterilizer chamber is smaller than that of the V-PRO maX 2 Sterilizer, therefore less sterilant is injected.
- V-PRO s2 Fast Cycle uses 1.1 g sterilant injection of VAPROX HC Sterilant (59% H₂O₂ solution) each sterilization pulse with 2 sterilization pulses per cycle, and a 50 °C chamber temperature. The volume of the V-PRO s2 sterilizer chamber is smaller than that of the V-PRO maX 2 sterilizer, therefore less sterilant is injected.
- The sterilizer physical components are the same between the V-PRO 1/1 Plus/maX/maX 2 Sterilizers except cabinetry differences for the V-PRO maX/maX 2 Sterilizer which have no impact on sterilization cycle execution or performance. The same chamber and sterilant delivery system is used for the V-PRO 1/1 Plus/maX/maX 2 Sterilizers. The V-PRO 60 chamber is smaller and the sterilant delivery system is the same as the V-PRO III Plus/maX/maX 2 Sterilizers.

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5.6 Controls after Sterilization Process

WARNING

Always wear protective gloves when emptying containers from the sterilizer.

Products should not undergo a change of environment with large temperature differences after sterilization process.

It should be expected to cool down to room temperature. Sweating will come to the fore in sudden temperature changes.

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Make sure the chemical indicator color on the paper labels has changed. (For paper labels see. 2.5)

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DE


Figure /Abbildung - 43

If a chemical-indicative safety seal is used prior to sterilization, the sterilization control can also be provided via the security seal.

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Figure /Abbildung - 44

In the paper filter models of sterilization containers, there are chemical indicators on the disposable filters. Internal control can be provided with disposable filters.

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Figure /Abbildung - 45

EN

The color conversion of the indicator color on the chemical indicator card after sterilization turns from Blue to Black.

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Figure /Figuur - 46

5.7 Transport and storage after sterilization

- Always carry the sterilization containers by their handles.
- Never carry sterilization container holding from the lid or lock system.
- Transport the sterilization container in such a way that mechanical damage will not occur.

Sterilization container storage room conditions ;

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Table / Tabelle -12

NOTE

The sterilization container systems have been validated as sterile for 6 Month and 1 Day under the appropriate storage conditions specified in table 12.

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DE

MALFUNCTION-A

Excessive amounts of condensate inside the sterilization container

Reason-A-1

Temperature of sterile materials too low prior to sterilization

Solution-A-1

Allow sterile materials to come to room temperature (appr. 20°C) before sterilization.

Reason-A-2

Textiles too damp

Solution-A-2

Sterilize dry textiles only

Reason-A-3

The sterilization container is overloaded

Solution-A-3

Load according to Table-9.

Reason-A-4

The containers to be sterilized are not correctly placed in sterilizer.

Solution-A-4

The heavier containers should be at the bottom and the total height should be 45 cm. Polypropylene containers should be preferred at the top

Reason-A-5

Improper sterilization program preferred.

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Solution-A-6

Observe the manufacturer's recommendations when using the autoclave. Check drying time and steam quality. If necessary, calibrate.

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Reason-A-7

The idle loop and vacuum test are not performed daily before starting sterilization.

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Solution-A-7

Perform idle cycle and vacuum tests before sterilization. Make sure your sterilizer is working properly.

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MALFUNCTION-B

Indicator color did not change after sterilization.

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Grund-B-1

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Solution-B-1

Have the autoclave repaired by the manufacturer.

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MALFUNCTION-C

Deformed the container after sterilization.

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Reason-C-1

Filter holes are closed during sterilization.

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Solution-C-1

Never close the holes in the body and the lid. Load according to 5.2

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Lid lock failure

Reason-D-1

Lock have been used instead of handles to carry the sterilization container system or the lids are not properly closed.

Solution-D-1

Please contact Manufacturer. Container with a broken lock is not suitable for sterilization.

Reason-D-2

The sterilization container system has had a physical impact on the lid locks.

Solution-D-2

Please contact Manufacturer. Container with a broken lock is not suitable for sterilization.

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WARNING NOTICE / WARNUNG

Before using first time

Our sterilization containers are generally delivered non-sterile. They must be cleaned, prepared and sterilized before being used for the first time!

Disposal

Please dispose of the sterilization container properly. Containers and accessories should be sterilized before disposal!

Incident

Any serious incident that has occurred in relation to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

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DE

7.0 Technical Service

For service and repairs, please contact REDA Instrumente GmbH

Modifications carried out on medical technical equipment may result in loss of warranty rights and forfeiture of applicable licenses.

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Web: www.REDA-INSTRUMENTE.DE

References Referenzen

- MDR 2017/745
- EN 285 in combination with EN ISO 17665
- EN ISO 11607
- EN 868-8
- DIN 58953-9
- DIN 58946-7
- HTM 01-01 Management and decontamination of surgical instruments used in acute care Part C and Part D
- CDC (Center for disease control and prevention)
- DAS (Disinfection Antiseptic Sterilization Association)

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Wenden Sie sich für Service und Reparaturen an REDA Instrumente GmbH